


Laryngeal Mask Disposable COLOUR CODED



22LMS10	SIZE 1	BLUE
22LMS15	SIZE 1.5	BLACK
22LMS20	SIZE 2	CLEAR
22LMS25	SIZE 2.5	GREEN
22LMS30	SIZE 3	ORANGE
22LMS40	SIZE 4	RED
22LMS50	SIZE 5	YELLOW

**REINFORCED
NOW AVAILABLE**

Laryngeal Mask Silicone Reusable

- 1) All silicone product with superior bio-compatibility, non reactive to patient tissue.
Flexible tube with a degree of hardness offers smooth insertion with maximum patient comfort.
- 2) Airway tube design allows safe use of suction catheters and bronchoscopes.
Reinforced rib design to inner tube makes for easier intubation.
- 3) Soft Silicone material on mask cushion is designed to perfectly conform with patient's hypopharynx.
- 4) Easy to read balloon and check valve identification. Also secure connection for excellent seal.
- 5) The device allows 40 uses after sterilisation.
- 6) Easy to use counting system "MEDICOUNT" allows simple scan of product to count its uses when it is in cleaning process. Once product reaches end of life the counting device indicator light flashes.

PRODUCT CODE NO	SIZE	I D (mm)	O D (mm)	MAXIMUM CUFF INFLATION VOLUME IN (MLS)	MAXIMUM ET TUBE SIZE (mm)	MAXIMUM FOB SIZE (mm)	PATIENT WEIGHT DETAILS
22LMS01RE	1	5.0	8.0	4	3.5	2.7	NEONATES / INFANTS UP TO 5KG
22LMS15RE	1.5	6.0	10.0	7	4.0	3.0	INFANTS 5 - 10KG
22LMS02RE	2	7.0	11.0	10	4.5	3.5	INFANTS / CHILDREN 10 - 20KG
22LMS25RE	2.5	8.5	13.0	14	5.0	4.0	CHILDREN 20 - 30KG
22LMS03RE	3	10.0	15.0	20	6.0 cuffed	5.0	CHILDREN 30 - 50KG
22LMS04RE	4	10.0	15.0	30	6.0 cuffed	5.0	ADULTS 50 - 70KG
22LMS05RE	5	11.0	16.5	40	7.0 cuffed	5.5	ADULTS 70 - 100KG

Laryngeal Mask Silicone Disposable

- 1) All silicone product with superior bio-compatibility, non reactive to patient tissue.
Flexible tube with a degree of hardness offers smooth insertion with maximum patient comfort.
- 2) Airway tube design allows safe use of suction catheters and bronchoscopes.
Reinforced rib design to inner tube makes for easier intubation.
- 3) Soft Silicone material on mask cushion is designed to perfectly conform with patient's hypopharynx.
- 4) Easy to read balloon and check valve identification. Also secure connection for excellent seal.
- 5) Colour coding to easily identify the balloon size being used.
- 6) Single Use Only

PRODUCT CODE NO	SIZE	I D (mm)	O D (mm)	MAXIMUM CUFF INFLATION VOLUME IN (MLS)	MAXIMUM ET TUBE SIZE (mm)	MAXIMUM FOB SIZE (mm)	PATIENT WEIGHT DETAILS
22LMS01	1	5.0	8.0	4	3.5	2.7	NEONATES / INFANTS UP TO 5KG
22LMS15	1.5	6.0	10.0	7	4.0	3.0	INFANTS 5 - 10KG
22LMS02	2	7.0	11.0	10	4.5	3.5	INFANTS / CHILDREN 10 - 20KG
22LMS25	2.5	8.5	13.0	14	5.0	4.0	CHILDREN 20 - 30KG
22LMS03	3	10.0	15.0	20	6.0 cuffed	5.0	CHILDREN 30 - 50KG
22LMS04	4	10.0	15.0	30	6.0 cuffed	5.0	ADULTS 50 - 70KG
22LMS05	5	11.0	16.5	40	7.0 cuffed	5.5	ADULTS 70 - 100KG

Laryngeal Mask Silicone Disposable REINFORCED

PRODUCT CODE NO	SIZE	I D (mm)	O D (mm)	MAXIMUM CUFF INFLATION VOLUME IN (MLS)	MAXIMUM ET TUBE SIZE (mm)	MAXIMUM FOB SIZE (mm)	PATIENT WEIGHT DETAILS
22LMR01	1	5.0	8.0	4	3.5	2.7	NEONATES / INFANTS UP TO 5KG
22LMR15	1.5	6.0	10.0	7	4.0	3.0	INFANTS 5 - 10KG
22LMR02	2	7.0	11.0	10	4.5	3.5	INFANTS / CHILDREN 10 - 20KG
22LMR25	2.5	8.5	13.0	14	5.0	4.0	CHILDREN 20 - 30KG
22LMR03	3	10.0	15.0	20	6.0 cuffed	5.0	CHILDREN 30 - 50KG
22LMR04	4	10.0	15.0	30	6.0 cuffed	5.0	ADULTS 50 - 70KG
22LMR05	5	11.0	16.5	40	7.0 cuffed	5.5	ADULTS 70 - 100KG