


Nasal Cannula


Oxygen flowrate 2 - 4 LPM - 24 - 30% O₂
 Oxygen flowrate 4 - 6 LPM - 30 - 40% O₂
1.8 Metre Oxygen Tubing included

Moulded Tips


Adult Flared Tips
221585




Adult Straight Tips
221582



Child Straight Tips
221601




Infant Straight Tips
221582




Softex Dipped Tips

Tips are dip moulded to give extra soft texture. Ideal for sensitive skin

Adult Curved Tips
221583



Adult Curved Tips
221601S



Adult Curved Tips
221602S

