



Medium Concentration Mask

Oxygen flowrate	2 - 4 LPM	-	24 - 30% O ₂
Oxygen flowrate	4 - 6 LPM	-	30 - 40% O ₂
Oxygen flowrate	6 - 8 LPM	-	40 - 50% O ₂

221420



Non-Rebreathing Mask

Hi-Concentration Oxygen delivery
with safety side vent.
Delivers up to 90% oxygen into mask
(Flowrate 8 - 15 LPM)

221423

221420 Adult Oxygen Mask

221470 Paediatric Oxygen Mask

221430 Adult Short Style Oxygen Mask

221480 Infant Oxygen Mask

221423 Adult Hi-Concentration Mask Non-Rebreathing

221425 Adult Hi-Concentration Mask Non-Rebreathing w/tubing