






GENERATOR		INSTRUMENT END CONNECTIONS				
5MM PIN CONNECTION >> Generators: Berchtold, Erbe, Martin & Eschmann 		cable length metres 3m SOCKET 48 10503048 10505048	SOCKET 40 10503040 10505040	SOCKET 30 10503030 10505030	SNARE 30 10503028 10505028	1.95MM PIN 10503095 10505095
3m AVAILABLE IN 3 & 5 METRE LENGTH CABLES 5m		FRONT & SIDE VIEWS OF INSTRUMENT END CONNECTIONS OLSEN 10503042 10505042		STORZ 10503010 10505010	STORZ 10503011 10505011	4MM PIN 10503044 10505044

S-CORD CABLES

GENERATOR	INSTRUMENT END CONNECTIONS			CLICK COUNTER
 S-CORD >> endoscopic earthing cables for valleylab generators for use in endoscopy	Olympus screw-fit for connecting short pin cable length metres 1m SOCKET 48 10401097	Same long pin SOCKET 30 10401098	Pentax screw-fit 4MM PIN 10401099	Accurate record of uses. 

ADAPTORS

GENERATOR ADAPTORS	
 8MM BOVIE CONNECTION >> with 4mm socket CODE >> 108040	CODE >> 108030  8mm Bovie connection with hook screw attachment

CODING FORMAT EXAMPLE

CABLE LENGTH: 3.0 METRES

1 0 4 0 3 0 4 8

MONOPOLAR CABLES | GENERATOR CONNECTION PIN SIZE | CABLE LENGTH | INSTRUMENT CONNECTION SOCKET STYLE

2 single use option please use prefix SU for single use

FOR ALL ORDERS CALL DIRECT **0700 - 4 - MEDITECH**