


# 12 >>

Cables bipolar

## GENERATOR




**TWIN 4MM PINS >>**  
compatible with most generators

## INSTRUMENT END CONNECTIONS

3m	EURO	STORZ	WOLF	US	AESCLAP (OLD)
5m	11443075 11445075	11443041 11445041	11443042 11445042	11443076 11445076	11443043 11445043

## GENERATOR




**COAXIAL PIN >>**  
Generators: Erbe, Downs & Siemens

## INSTRUMENT END CONNECTIONS

3m	EURO	STORZ	WOLF	US	AESCLAP (OLD)
5m	11423075 11425075	11423041 11425041	11423042 11425042	11423076 11425076	11423043 11425043

## GENERATOR



**COAXIAL PIN >>**  
Generators: Berchtold, Wolf, & Martin

## INSTRUMENT END CONNECTIONS

3m	EURO	STORZ	WOLF	US	AESCLAP (OLD)
5m	11453075 11455075	11453041 11455041	11453042 11455042	11453076 11455076	11453043 11455043

## GENERATOR



**LEMO CONNECTOR >>**  
Generators: Older Eschmann (Iemo)

## INSTRUMENT END CONNECTIONS

3m	EURO	STORZ	WOLF	US	AESCLAP (OLD)
5m	11433075 11435075	11433041 11435041	11433042 11435042	11433076 11435076	11433043 11435043

FOR ALL ORDERS CALL DIRECT **0700 - 4 - MEDITECH**

TEL No +44 (0) 1258 471 770  
FAX No +44 (0) 1258 471 772

Meditech Systems Limited  
Unit 3 Richmar  
Butt's Pond Ind.Estate  
Sturminster Newton  
Dorset DT10 1AZ

