



16 >>

Disposable bipolar forceps

SINGLE USE FORCEPS

INTEGRAL CABLE

FORCEPS LEG



TIP SIZES COLOUR CODED

FEATURES:

- STAINLESS STEEL LEGS
 - Providing same feel as reusable
- NYLON COATING
 - Same insulation as reusable
- INTEGRAL CABLE
 - Cost effective, choice of pins

SINGLE USE BIPOLAR FORCEPS

STYLE	LENGTH	TIP SIZES		
		2MM	1MM	0.5MM
STRAIGHT				
Cawthorne	80mm	su762008044	su761008044	su760508044*
Cawthorne	110mm	su762011044	su761011044	su760511044
Standard pattern	110mm	su772011044	su771011044	su770511044
Standard pattern	130mm	su772013044	su771013044	su770513044
Standard pattern	150mm	su772015044	su771015044	su770515044
Standard pattern	180mm	su772018044	su771018044	su770518044
Standard pattern	200mm	su772020044	su771020044	su770520044
Adson pattern	130mm	su792013044	su791013044	su790513044
Adson pattern	150mm	su792015044	su791015044	su790515044


*PLEASE NOTE:

for Kocher tooth tips add the letter "K" (eg 762008044K)
for serrated tips add the letter "S" (eg 762008044S)

STYLE	LENGTH	TIP SIZES		
		2MM	1MM	0.5MM
ANGLED				
Bayonet	150mm	su782015044	su781015044	su780515044
Bayonet	200mm	su782020044	su781020044	su780520044
Angled tip	80mm	su752008044	su751008044	su750508044
Angled tip	110mm	su752011044	su751011044	su750511044
Angled tip	130mm	su752013044	su751013044	su750513044
Angled tip	150mm	su752015044	su751015044	su750515044
Angled tip	180mm	su752018044	su751018044	su750518044
Angled tip	200mm	su752020044	su751020044	su750520044

*PLEASE NOTE:

for Kocher tooth tips add the letter "K" (eg 762008044K)
for serrated tips add the letter "S" (eg 762008044S)



+44 (0) 1258 471 770
+44 (0) 1258 471 772

TEL: 01258 471 770
FAX: 01258 471 772

Meditech Systems Limited
Unit 3 Richmar
Butt's Pond Ind.Estate
Sturminster Newton
Dorset DT10 1AZ

meditech
systems limited

FOR ALL ORDERS CALL DIRECT 0700 - 4 - MEDITECH