

"INSU-TIP" FORCEPS


FEATURES:

INSULATED TIP

Outer surfaces of tips are insulated to prevent burning to peripheral tissue

INSULATED TIP



70% reduction of exposed metal on tips.



TIP SIZES COLOUR CODED

INSULATED TIP

Energy focused between tips allowing lower power setting to be used



"INSU-TIP" BIPOLAR FORCEPS				
STYLE	LENGTH	TIP SIZES		
		2MM	1MM	0.5MM
STRAIGHT				
Cawthorne	80mm	su862008044	su861008044	su860508044*
Cawthorne	110mm	su862011044	su861011044	su860511044
Standard pattern	110mm	su872011044	su871011044	su870511044
Standard pattern	130mm	su872013044	su871013044	su870513044
Standard pattern	150mm	su872015044	su871015044	su870515044
Standard pattern	180mm	su872018044	su871018044	su870518044
Standard pattern	200mm	su872020044	su871020044	su870520044
Adson pattern	130mm	su892013044	su891013044	su890513044
Adson pattern	150mm	su892015044	su891015044	su890515044
*PLEASE NOTE: for Kocher tooth tips add the letter "K" (eg 762008044K) for serrated tips add the letter "S" (eg 762008044S)				
STYLE	LENGTH	TIP SIZES		
		2MM	1MM	0.5MM
ANGLED				
Bayonet	150mm	su882015044	su881015044	su880515044
Bayonet	200mm	su882020044	su881020044	su880520044
Angled tip	80mm	su852008044	su851008044	su850508044
Angled tip	110mm	su852011044	su851011044	su850511044
Angled tip	130mm	su852013044	su851013044	su850513044
Angled tip	150mm	su852015044	su851015044	su850515044
Angled tip	180mm	su852018044	su851018044	su850518044
Angled tip	200mm	su852020044	su851020044	su850520044
*PLEASE NOTE: for Kocher tooth tips add the letter "K" (eg 762008044K) for serrated tips add the letter "S" (eg 762008044S)				