


























30 >>
Light cables

INSTRUMENT CONNECTION


AVAILABLE WITH FIXED OR REMOVABLE END CONNECTORS

LIGHT CABLES						
LIGHTSOURCE CONNECTION		INSTRUMENT CONNECTION		DESCRIPTION		
	1		1	LIGHTSOURCE CONNECTORS 1. ACMI (AFI, Aspen, Cuda, Eder, KLI, Lutex, Pilling, Stryker, V.mueller, Wehmer, Z immer). 2. WOLF DYONICS 3. WOLF EXTENDED 4. KARL STORZ OLYMPUS (5.0 MAX. BUNDLE DIA) 5. OLYMPUS 6. PILLING - 6 (3.5 MAX. BUNDLE DIA) 7. DESIGNS FOR VISION 8. MACHIDA 9. PENTAX INSTRUMENT CONNECTORS 1. ACMI (AFI, Aspen, Cuda, Eder, KLI, Lutex, Pilling, Stryker, V.mueller, Wehmer, Z immer). 2. WOLF DYONICS 3. KARL STORZ OLYMPUS ACMI-1 (AFI, Aspen, Cuda, Eder, KLI, Lutex, Pilling, Stryker, V.mueller, Wehmer, Z immer). 4. WOLF DYONICS CLIP 5. SCHOLLY 6. KARL STORZ OLYMPUS (5.0 MAX. BUNDLE DIA) 7. PILLING 8. DESIGNS FOR VISION 9. ACMI MALE 10. ACMI CLIP-ON 11. ACMI 45° MALE 12. ACMI 90° MALE		
	2		2			
	3		3			
	4		4			
	5		5			
	6		6			
	7		7			
	8		8			
	9		9			
	10		10			
	11		11			
	12		12			

FIBREOPTIC BUNDLE SIZE (mm)





CABLE END VIEW



Available from 2mm to 6.5mm diameter

CABLE LENGTHS: 180/225/275CM

CODING FORMAT EXAMPLE



+44 (0) 1258 471 770
+44 (0) 1258 471 772

Meditech Systems Limited
Unit 3 Richmar
Butt's Pond Ind.Estate
Sturminster Newton
Dorset DT10 1AZ



FOR ALL ORDERS CALL DIRECT **0700 - 4 - MEDITECH**