



**LIGHT CARRIERS**




LIGHT CARRIER FOR RHYSEVANS MICRO LARYNGOSCOPE  
**CODE >> 88C7500**


Example of light carrier fitting to Rhys - Evans micro laryngoscope  
Custom styles available



LIGHT CARRIER FOR USE WITH NEGUS LARYNGOSCOPE  
**CODE >> 88C700**



LIGHT CARRIER FOR USE WITH NEGUS OESOPHAGOSCOPES  
**CODE >> 88C1000**




LIGHT CARRIER FOR USE WITH RHYSEVANS MICRO LARYNGOSCOPE  
**CODE >> 88C7500**


FIBRE LOADING REPAIR SERVICE



**SCOPE REPAIR SERVICE**




LIGHT CARRIER FOR USE WITH RNTNE BLACKENED LARYNGOSCOPES  
**CODE >> 88C4300**




LIGHT CARRIER FOR USE WITH NEGUS BRC  
**CODES:**

- 88C100** 40cm (2mm dia)
- 88C200** 30cm (2mm dia)
- 88C300** 27.5cm (2mm dia)
- 88C400** 45cm (2mm dia)
- 88C500** 27.5cm (1.6mm)
- 88C600** 18cm (2mm dia)



**LIGHT CARRIERS**

CODE	SIZE	COMPATIBLE INSTRUMENT
fibre light carriers for bronchoscopes		
<b>88C100</b>	40cm	SS4085-SS3042-SS4024-SS4281-SS4025*
<b>88C200</b>	30cm	SS4026-SS4027
<b>88C300</b>	27.5cm	SS4028
<b>88C400</b>	45cm	SS4029
<b>88C500</b>	27.5cm	SS4030-SS4280
<b>88C600</b>	18cm	SS4670
Twin fibre light carriers for laryngoscopes		
<b>88C700</b>	—	SS4083-SS4021-SS4023*
<b>88C7A00</b>	—	SS4084
<b>88C800</b>	—	SS4020-SS4022
<b>88C8A00</b>	—	SS4102
<b>88C900</b>	—	SS4101
<b>88C9A00</b>	—	SS4103
<b>88C1000</b>	—	SS4086-SS4031*
<b>88C10A00</b>	—	SS4106
<b>88C10B00</b>	—	SS4032
<b>88C1100</b>	—	SS4520-SS4521
<b>88C11A00</b>	—	SS4033-SS4522
<b>88C1200</b>	—	SS4523
<b>88C1300</b>	—	SS4524



**LIGHT CARRIERS**

Example of light carrier fitting to proctoscope  
Custom styles available

**CODE >> 88C1600**

\*PLEASE NOTE:  
Above fibre light carriers will fit both polished metal and stainless steel instruments